




## Pycnometer - Cells/Spheres/Sleeves



### **LARGE 135cm<sup>3</sup>**

SAMPLE CELL 04000-0414  
SPHERE p/n 01500-LARGE  
SPHERE VOLUME 56.5592cc  
Dia 47.6mm / 1.874in

### **MEDIUM 50cm<sup>3</sup>**

SAMPLE CELL 04000-3314  
SPHERE p/n 01500-MEDIUM  
SPHERE VOLUME 28.9583cc  
Dia 38.07mm / 1.499in  
MEDIUM SLEEVE 04000-3316

### **SMALL 10cm<sup>3</sup>**

SAMPLE CELL 04000-3319  
SPHERE p/n 01500-SMALL  
SPHERE VOLUME 7.0699cc  
Dia 23.79mm / 0.936in  
SMALL SLEEVE 04000-3318

### **MICRO 4.5cm<sup>3</sup>**

MICRO CELL 04000-3320  
SPHERE p/n 01500-MICRO  
SPHERE VOLUME 1.0725cc  
Dia 12.68MM / 0.5in  
MICRO SLEEVE 04000-3321

### **MESO 1.8 cm<sup>3</sup>**

MESO CELL 04000-5418  
SPHERE p/n 01500-MICRO  
SPHERE VOLUME 1.0725cc  
Dia 12.68MM / 0.5in  
MICRO SLEEVE 04000-3321

### **Nano 0.5cm<sup>3</sup>**

Nano CELL 04000-3707  
SPHERE p/n 01500-NANO  
SPHERE VOLUME 0.09cc  
Dia 5.55MM / 0.218 in  
MICRO SLEEVE 04000-3321  
MICRO CELL 04000-3320

CELL SPHERES SLEEVES